VITALPEAK MSG 150

ENG INSTRUCTIONS FOR USE

VITALPEAK MSG 150

Back and neck massage device

Thank you for purchasing the VITALPEAK product

Safety when working with the device and use. Read carefully the enclosed safety instructions.

- For proper handling of the product, read the latest version of the manual, which can be downloaded at: www.vitalpeak.info
- The device is intended for use in a domestic or private environment only, not for commercial or medical purposes
 - You can indulge in a pleasant massage with this massage device without outside help or use the device during / after exercise or other physical exertion.
 - The manufacturer is not liable for damage caused by improper or incorrect use.

Using the device:

1. Connect the included adapter to the device
2. Plug the adapter into the electrical outlet
3. Briefly press the power button to start the massage
4. If you also want to switch on the heating function, briefly press the heating button
5. If you want to change the direction of the massage, press the change direction button
6. Briefly press the power button to turn off
7. Unplug the adapter after turning off

 Heating button

 Power button

 Change direction button

Automatic shutoff

The massage should not last more than 20 minutes. The device shuts off automatically after 20 minutes. Wait at least 20 minutes before next massage. Too long a massage can lead to stretching and excessive stimulation of muscles and nerves.

Specifications:

Adapter input specifications: 100 – 240 V AC
Frequency: 50 – 60 Hz
Adapter output specifications: 13 V DC, 1.8 A
Power: 24 W
Automatic shutoff: 20 minutes

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